RIGENERA

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Sonja Ellen Lobo and Treena Livingston Arinze. Biphasic Calcium Phosphate Ceramics for Bone Regeneration and Tissue Engineering Applications. *Review. Materials* 2010, 3, 815-826;

Gauthier, O.; Bouler, J.M.; Aguado, E.; LeGeros, R.Z.; Pilet, P.; Daculsi, G. Elaboration conditions influence physicochemical properties and in vivo bioactivity of macroporous biphasic calcium phosphate ceramics. *J. Mater. Sci. Mater. Med.* 1999, 10, 199-204.

Barrère, F.; van Blitterswijk, C.A.; de Groot, K. Bone regeneration: Molecular and cellular interactions with calcium phosphate ceramics. *Int. J. Nanomed.* 2006, 1, 317-332

Le Nihouannen, D.; Duval, L.; Lecomte, A.; Julien, M.; Guicheux, J.; Daculsi, G.; Layrolle, P. Interactions of total bone marrow cells with increasing quantities of macroporous calcium phosphate ceramic granules. *J. Mater. Sci. Mater. Med.* 2007, 18, 1983-1990.

Khan, Y.; Yaszemski, M.J.; Mikos, A.G.; Laurencin, C.T. Tissue engineering of bone: Material and matrix considerations. *J. Bone Joint. Surg. Am.* 2008, 90, 36-42.

Davies, J.E. *Bone Engineering*; Em Squared Inc: Toronto, Canada, 2000; pp. 516-525.

Hing, K.A.; Wilson, L.F.; Buckland, T. Comparative performance of three ceramic bone graft substitutes. *Spine J.* 2007, 7, 475-490.

Adamopoulos, O.; Papadopoulos, T. Nanostructured bioceramics for maxillofacial applications. *J. Mater. Sci. Mater. Med.* 2007, 18, 1587-1597.

LeGeros, R.Z.; Lin, S.; Rohanzadeh, R.; Mijares, D.; LeGeros, J.P. Biphasic calcium phosphate bioceramics: Preparation, properties and applications. *J. Mater. Sci. Mater. Med.* 2003, 14, 201-209.

Daculsi, G.; LeGeros, R. Biphasic calcium phosphate (BCP) bioceramics: Chemical, physical and biological properties. *Enc. Biomat. Biomed. Eng.* 2006, 1-9

Habibovic, P.; Sees, T.M.; van den Doel, M.A.; van Blitterswijk, C.A.; de Groot, K. Osteoinduction by biomaterials—physicochemical and structural influences. *J. Biomed. Mater. Res. A* 2006, 77, 747-762.

Lobo, S.E.; Wykrota, F.H.; Oliveira, A.C.; Kerkis, I.; Mahecha, G.B.; Alves, H.J. Quantification of bone mass gain in response to the application of biphasic bioceramics and platelet concentrate in critical-size bone defects. *J. Mater. Sci. Mater. Med.* 2009, 20, 1137-1147.

Daculsi, G.; LeGeros, R.Z.; Nery, E.; Lynch, K.; Kerbel, B. Transformation of biphasic calcium phosphate ceramics in vivo: Ultrastructural and physicochemical characterization. *J. Biomed. Mater. Res.* 1989, 23, 883-894.

Daculsi, G.; LeGeros, R.Z.; Heughebaert, M.; Barbieux, I. Formation of carbonate-apatite crystals after implantation of calcium phosphate ceramics. *Calcif. Tissue Int.* 1990, 46, 20-27.

ASTM F 1185, Standard Specification for Composition of Hydroxyapatite for Surgical Implants

ASTM F 1088, Standard Specification for Beta-Tricalcium Phosphate for Surgical Implantation

Daculsi G. et al. Spongy and Cortical bone substitution kinetics at the expense of the macroporous biphasic calcium phosphate: animal and human evidence. *Bioceramics*. 1999, 12 pp. 287-290

Ducheyne P. et al. The effect of calcium phosphate ceramic composition and structure on in vivo behaviour. I. Dissolution. *J. Biomed. Mater. Res.* 1993, 27 pp. 25-34

Daculsi G, Jegoux F, Layrolle P. The micro macroporous biphasic calcium phosphate concept for bone reconstruction and tissue engineering. In: Basu B, Katti DS, Kumar A, editors. *Advance biomaterials: fundamentals processing and applications*. New York: Wiley-American Ceramic Society; 2009. p. 101-42

Daculsi G, Baroth S, LeGeros RZ. 20 years of biphasic calcium phosphate bioceramics development and applications. In: Narayan R, Colombo P, Singh D, Salem J, editors. *Advance in bioceramics and porous ceramics II*. New York: Wiley-American Ceramic Society; 2010. p. 45-58.

Daculsi G. Biphasic calcium phosphate concept applied to artificial bone, implant coating and injectable bone substitute. *Biomaterials* 1998;19: 1473-8.

Lecomte A, Gauthier H, Bouler JM, Gouyette A, Pegon Y, Daculsi G, et al. Biphasic calcium phosphate: a comparative study of interconnected porosity in two ceramics. *J Biomed Mater Res B Appl Biomater* 2008;84B:1-6.

Tancret F, Bouler JM, Chamousset J, Minois LM. Modelling the mechanical properties of microporous and macroporous biphasic calcium phosphate bioceramics. *J Eur Ceram Soc* 2006;26:3647-56.

Daculsi G, Weiss P, Bouler JM, Gauthier O, Millot F, Aguado E. Biphasic calcium phosphate/hydroxylpolymer composites: a new concept for bone and dental substitution biomaterials. *Bone* 1999;25(Suppl. 2):595-615.

Gauthier O, Goyenvalle E, Bouler JM, Guicheux J, Pilet P, Weiss P, et al. Macroporous biphasic calcium phosphate ceramics versus injectable bone substitute: a comparative study 3 and 8 weeks after implantation in rabbit bone. *J Mater Sci: Mater Med* 2001;12:385-90.

Daculsi G. Biphasic calcium phosphate granules concept for injectable and mouldable bone substitute. *Adv Sci Technol* 2006;49:9-13.

Mangano C, Sinjari B, Shibli JA, Mangano F, Hamisch S, Piattelli A, Perrotti V, DDS, Giovanna Iezzi G; A Human Clinical, Histological, Histomorphometrical, and Radiographical Study on Biphasic HA-Beta-TCP 30/70 in Maxillary Sinus Augmentation; *Clinical Implant Dentistry and Related Research*, Volume 17, Number 3, 2015.

Uzeda M, J, Rodrigo Figueiredo de Brito Resende, Suelen Cristina Sartoretto, Adriana Terezinha Neves Novellino Alves, Jose Mauro Granjeiro, Monica Duanas Calasans-Maia; Randomized clinical trial for the biological evaluation of two nanostructured biphasic calcium phosphate biomaterials as a bone substitute; *Clin Implant Dent Relat Res*. 2017;19:802-811

Corbella S, Taschieri S, Weinstein R, Del Fabbro M.; "Histomorphometric outcomes after lateral sinus floor elevation procedure: a systematic review of the literature and meta-analysis". *Clin. Oral Impl. Res.* 00, 2015, 1-17 doi: 10.1111/dr.12702

W. L. Grayson, M. Fröhlich, K. Yeager et al., "Engineering anatomically shaped human bone grafts," *Proceedings of the National Academy of Sciences of the United States of America*, vol. 107, no. 8, pp. 3299-3304, 2010.

J. Li, L. Zhang, S. Lv, S. Li, N. Wang, and Z. Zhang, "Fabrication of individual scaffolds based on a patient-specific alveolar bone defect model," *Journal of Biotechnology*, vol. 151, no. 1, pp. 87-93, 2011.

J. Ruestmeyer, A. Busch, and A. Sari-Rieger, "Application of computer-aided designed/computer-aided manufactured techniques in reconstructing maxillofacial bony structures. *Oral and Maxillofacial Surgery*, vol. 18, no. 4, pp. 471-476, 2014.

F. Mangano, P. Zecca, S. Pozzi-Taubert et al., "Maxillary sinus augmentation using computer-aided design/computer-aided manufacturing (CAD/CAM) technology," *The International Journal of Medical Robotics and Computer Assisted Surgery*, vol. 9, no. 3, pp. 331-338, 2013.

R. Kontio, "Update on mandibular reconstruction: computer aided design, imaging, stem cells and future applications," *Current Opinion in Otolaryngology & Head and Neck Surgery*, vol. 22, no. 4, pp. 307-315, 2014.

M. Esposito, M. G. Grusovin, P. Felice, G. Karatzopoulos, H.V. Worthington, and P. Coulthard, "The efficacy of horizontal and vertical bone augmentation procedures for dental implants—a Cochrane systematic review," *European Journal of Oral Implantology*, vol. 2, no. 3, pp. 167-184, 2009.

B. Al-Nawas and E. Schiegnitz, "Augmentation procedures using bone substitute materials or autogenous bone—a systematic review and meta-analysis," *European Journal of Oral Implantology*, vol. 7, supplement 2, pp. s219-s234, 2014.

Tarsitano, S. Mazzoni, R. Cipriani, R. Scotti, C. Marchetti, and L. Ciocca, "The CAD-CAM technique for mandibular reconstruction: an 18 patients oncological case-series," *Journal of Cranio-Maxillofacial Surgery*, vol. 42, no. 7, pp. 1460-1464, 2014.

Y.-F. Liu, L.-W. Xu, H.-Y. Zhu and S. S.-Y. Liu, "Technical procedures for template-guided surgery for mandibular reconstruction based on digital design and manufacturing," *Biomedical Engineering Online*, vol. 13, no. 63, pp. 1-15, 2014.

F. Tamimi, J. Torres, K. Al-Abedalla et al., "Osseointegration of dental implants in 3D-printed synthetic onlay grafts customized according to bone metabolic activity in recipient site," *Biomaterials*, vol. 35, no. 21, pp. 5436-5445, 2014.

F. Mangano, A. Macchi, J.A. Shibli et al., "Maxillary ridge augmentation with custom-made CAD/CAM scaffolds. A 1-year prospective study on 10 patients," *Journal of Oral Implantology*, vol. 40, no. 5, pp. 561-569, 2014.

M. Figliuzzi, F. G. Mangano, L. Fortunato et al., "Vertical ridge augmentation of the atrophic posterior mandible with custom made, computer-aided design/computer-aided manufacturing porous hydroxyapatite scaffolds," *Journal of Craniofacial Surgery*, vol. 24, no. 3, pp. 856-859, 2013.

Ella, M. Laurentjoye, C. Sedarat, J.-C. Coutant, E. Masson, and A. Rouas, "Mandibular ridge expansion using a horizontal bone-splitting technique and synthetic bone substitute: an alternative to bone block grafting?" *The International Journal of Oral & Maxillofacial Implants*, vol. 29, no. 1, pp. 135-140, 2014.

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